

# Trail Handles Installation Guide

---

Have a question? Call 877-345-4748 or email [support@gogrit.us](mailto:support@gogrit.us).

## Installation


Loosen the nylon of the seat back. Then, separate the two halves of each Trail Handle Clamp by removing the four front-facing screws. Noting left and right clamp setup, attach the Trail Handle Clamps **2.8"** from the tops of the canes (this positions the clamps between two nylon straps). Reaffix the seat back nylon.

**Note:** Loosening all four straps will allow the top set to reach over the Trail Handle Clamps.




## Adjustments

A) To straighten the Trail Handles or move them up or down the canes, loosen the four front-facing clamp bolts. Once loosened, slide or adjust each handle to your desired position and retighten the four bolts.



B) To adjust the angle of the Trail Handles, loosen the two pinch bolts on the rear-facing half of the Trail Handle Clamps. Rotate



## Note:

When tightening the four Trail Handle Clamp bolts, it is best to use a star pattern, as though tightening the lug nuts of a car.

If either Trail Handle rotates during use, apply slightly more torque, evenly, to each of the four Trail Handle Clamp bolts and/or the two pinch bolts, depending on the direction of the rotation.

